

Experiment No.7

To estimate Copper in given sample Calorimetrically/
Spectrometrically




Experiment No.7 To estimate Copper in given sample Calorimetrically

- **AIM :** To estimate Copper in given sample Calorimetrically
- **APPARATUS :** Colorimeter/ Spectrometer, volumetric flasks etc.
- **REAGENTS:** Copper sulphate solution, Ammonia solution



Experiment No.7

To estimate Copper in given sample Calorimetrically




- The working of colorimeter/ Spectrometrically is based on Beer's- Lambert's law.
- It state that,
- When beam of monochromatic light passing throw homogeneous absorbing medium, decrease in intensity of light is directly proportional to intensity of incident light as well as concentration of absorbing medium



Experiment No.7

To estimate Copper in given sample Calorimetrically




- **Transmission Or Transmittance: (T):** It is simple ratio of intensity of transmitted light to incident light
i.e. $T = I_t / I_0$
- **Optical density (OD) or Absorbance (A):** It is simple logarithmic ratio of intensity of incident light to transmitted light. $OD \text{ or } A = \log I_0 / I_t$



<https://youtu.be/wxrAELeXlek> link for UV-visible Spectrophotometer

Experiment No.7

To estimate Copper in given sample Calorimetrically



Experiment No.7

To estimate Copper in given sample Calorimetrically

- **Applications of Colorimetry**
- Concentration of solution as low as 10^{-7} M which cannot be determined by the conventional method like volumetric & gravimetric.
- Determination of iron.
- Determination of ammonia using Nessler's reagent.
- Determination of aluminum using Erio-chrome cyanine
- Determination of arsenic by molybdenum blue method.
- Simultaneous determination of chromium and manganese in steel.
- Analysis of ore, minerals, alloys & other industrial raw materials & finished products.
- Analysis of environmental samples.
- In determination of molecular composition of complexes.
- For studying Cis & Trans isomers.



Experiment No.7

To estimate Copper in given sample Calorimetrically

- **PROCEDURE:**

- **5.** Make a blank solution by diluting 5ml NH_3 solution in a 50 ml volumetric flask. Add distilled water up to the mark & mix well.
- **6.** Selection of λ_{max} or filter: **a)** Select filter no.1, Insert blank solution test tube in a sample holder & Press % T button. Followed by pressing unlock button & then calibrate button. Observe 100 on display. Immediately remove blank test tube and insert test tube containing 100ppm solution from flask 5. Press **ABS** button and record reading in table 1. Repeat above procedure for all filters, and record readings.
- From table no. 1 find out filter at which **ABS** is maximum. Select that filter as λ_{max} .



Experiment No.7

To estimate Copper in given sample Calorimetrically

- **PROCEDURE:**

- **7.** Select this filter, calibrate with blank solution and record **ABS** for all solutions.(2, 3, 4, 5, and unknown) using test tube used for blank solution. Record readings in table 2. .
- **c)** Plot a graph between ABS (OD) & Concentration. Find out unknown concentration graphically.
 - X axis: Concentration of Cu in mg Y axis: OD
 - Nature of graph: Straight line passing through origin.



Experiment No.7

To estimate Copper in given sample Calorimetrically

OBSERVATION TABLE

Table No.1: For λ max

Sr. No.	Filter (nm)	Optical Density (ABS)
1	400	0.01
2	420	0.01
3	470	0.04
4	500	0.08
5	530	0.08
6	620 λ max	0.14
7	660	0.13
8	700	0.06



Experiment No.7

To estimate Copper in given sample Calorimetrically

OBSERVATION TABLE

Table2 : For OD at 620 λ max

Vol. of std. Sol	Conc. mg/l	OD at 620 λ max
05ml	25	0.04
10ml	50	0.08
15ml	75	0.12
20ml	100	0.16
Test Solution	Unknown	0.10
Conc. Of Cu in unknown solution from Graph		
62.5 mg/l		



Experiment No.7
To estimate Copper in given sample Calorimetrically

Conc. mg/l	OD at 620 λ max
0	0
25	0.04
50	0.08
75	0.12
100	0.16
Unknown	0.1

