

pH Meter Setup



EXPERIMENT NO 6

To perform acid base titration pH metrically

- **AIM:** To perform Strong acid Vs Strong base titration pH metrically
- **APPARATUS:** Glass electrode, beaker, pH meter etc.
- **REAGENTS:** Given HCl acid (Bottle A) and 0.1N NaOH (Bottle B)



EXPERIMENT NO 6


To perform acid base titration pH metrically


- **THEORY:** The modern pH meter is an electronic digital voltmeter scaled to read pH directly and may range from a comparatively simple handy instrument suitable for use in the field also. It is useful technique to find out
 - 1. Equivalent point of any type of Acid – Base titrations
 - 2. To carry out above titrations in presence of strong oxidant, reductants, viscous protein mediums, as well as colored, colloidal, complex, poisonous solutions.
 - 3. To carry out weak acid weak base titrations or for dilute solutions.
 - 4. To carry out equivalent point of mixtures of two or more acids or bases in single titration.
 - 5. To find out equivalence point of di & tri basic or acidic substances in single titrations.
 - 6. To find out equivalence point of acid base titration in non aqueous medium.

EXPERIMENT NO 6

To perform acid base titration pH metrically


SCHEMATIC DIAGRAM OF pH METER





EXPERIMENT NO 6

To perform acid base titration pH metrically



- **Glass electrode:** Lower bulb of glass electrode is highly sensitive to H^+ ion concentration. It is made up of Corning glass. Thickness of the walls of glass bulb varies from 0.001 – 0.05 mm. It possesses low M.P. & high electrical conductivity. It is filled with 0.1M HCl solution and Ag wire coated with AgCl immerse in it. It is represented as Pt or Ag/ AgCl / 0.1 N HCl / Glass bulb/ Test solution.
- When glass electrode immersed in unknown solution, due to variation in concentration of H^+ ions inside & out side of the bulb potential is developed on inner & outer wall of glass bulb due to exchange of H^+ & free Na^+ present in crystal lattice of glass bulb, denoted by E1 & E2 potential. Difference between this two potential is called boundary potential. This potential difference is measured with potentiometer. This is displayed as pH on display of digital pH meter

EXPERIMENT NO 6

To perform acid base titration pH metrically

- **Procedure for pH- metric titration:-**

- 1) Clean the electrode by distilled water & Calibrate equipment.
- 2) Take 25 ml of given acid solution by pipette in a 50 ml beaker and immerse the electrode. Note down the pH meter reading for zero ml base addition.
- 3) Fill the burette with standard base solution (0.1N NaOH). Remove air bubble and adjust zero level.
- 4) Add Standard base solution gradually from the burette i.e. 0.5 ml at a time to acid solution, stir the same for 1 minute. Note the corresponding value of pH. Record readings in table given


EXPERIMENT NO 6: To perform acid base titration pH metrically

Volume of 0.1N NaOH added in ml	Observed pH	Δv	ΔpH	$\Delta pH / \Delta v$	$\Delta^2 pH / \Delta v^2$
0	1.87	---	-----	-----	-----
1	1.88	1.0 (1.0-1)	0.01 (1.88-1.87=0.01)	0.01 (0.01/1.0)	0.04 (0.05-0.01=0.04)
2	1.93	1.0	0.05 (1.93-1.88=0.05)	0.05	0.01
3	1.99	1.0	0.06	0.06	0.05
4	2.10	1.0	0.11	0.11	0.09
5	2.31	1.0	0.20	0.20	0.33
6	2.88	1.0	0.53	0.53	7.57
7	10.98	1.0	8.10	8.10	-7.46
8	11.62	1.0	0.64	0.64	-0.42
9	11.84	1.0	0.22	0.22	-0.08
10	11.98	1.0	0.14	0.14	

EXPERIMENT NO 6

To perform acid base titration pH metrically

Volume of 0.1N NaOH added in ml	Observed pH
0	1.87
1	1.88
2	1.93
3	1.99
4	2.1
5	2.31
6	2.88
7	10.98
8	11.62
9	11.84
10	11.98



8

9

10

11

12

13

EXPERIMENT NO 6

To perform acid base titration pH metrically

Volume of 0.1N NaOH added in ml	$\Delta pH / \Delta v$
0	0
1	0.01
2	0.05
3	0.06
4	0.11
5	0.2
6	0.53
7	8.1
8	0.64
9	0.22
10	0.14

$\Delta pH / \Delta v$ Vs ml of 0.1N NaOH added

Eq. point 'b' = 7.0 ml

EXPERIMENT NO 6

To perform acid base titration pH metrically

CALCULATION:

Equivalence point from graph. a) = 6.5 ml b) = 7.0 ml c) = 6.5 ml

Normality of Acid = $N_1V_1 = N_2V_2 = 0.1 \times 6.5 / 25$

= 0.026 N

Result: 1) Normality of Acid 0.026 N

2) pH range at equivalence point 2.88 to 10.98

3) Type of acid base titration Strong Acid Vs Strong Base