

Presenting...



Give control



Stop presenting

Experiment No 8

**To determine viscosity of given liquid by
Ostwald's Viscometer**



Type here to search



Presenting...



Give control



Stop presenting

Experiment No 8

To determine viscosity of given liquid by Ostwald's Viscometer

- **AIM:** To determine viscosity of given liquid.
- **APPARATUS:** Ostwald's Viscometer, Stop watch, Air dryer etc
- **REAGENTS:** Given liquid, Distilled water, Acetone etc



Type here to search



Presenting...



Give control



Stop presenting

Experiment No 8

To determine viscosity of given liquid by Ostwald's Viscometer

• THEORY:

- It is property of liquid or fluid by virtue of which it offers resistance to its own flow.
- The flowing liquid consists of a series of parallel layers moving one over other.
- Any two layers will move with velocity of different magnitudes.
- Top layer moves faster than the next lower layer, due to internal friction called **viscous drag**.
- The unit of viscosity is poise & is define as, when force of one dyne is required to maintain a relative velocity difference of 1cm/sec between two parallel layers, separated by 1cm
- The coefficient of viscosity of a liquid is 1 poise. A smaller corresponding unit is centipoise.
- **1poise = 100 centipoise.**




Type here to search



Experiment No 8

To determine viscosity of given liquid by Ostwald's Viscometer



- Simple method for estimate viscosity of any liquid at room temperature is **Ostwald's Viscometer**.
- In this method fixed volume of liquid & pure water allow to flow through a standard capillary tube in a fixed distance. The time in second taken by liquid & pure water is a comparative parameter to find out viscosity of unknown liquid at room temperature



Type here to search



Experiment No 8

To determine viscosity of given liquid by Ostwald's Viscometer

OBSERVATION TABLE

Name of Liquid	Sp. Gravity	Time required in sec.	Mean time in sec
Distilled water	0.997	1 21.94 2 21.37 3 21.13	21.48
Liquid Turpentine	0.852	1 26.33 2 25.91 3 25.75	25.99



Experiment No 8

To determine viscosity of given liquid by Ostwald's Viscometer

CALCULATIONS

$$\text{Viscosity of liquid } \eta = \frac{d_o * t_o}{d_w * t_w} \times \eta_w \text{ Centipoise}$$

Where, η_w = Viscosity of Water (0.8937), d_o = density of Turpentine, d_w = density of water, t_o = time
For Turpentine, t_w = time for water.

$$\text{Viscosity of liquid } \eta = (0.852 \times 25.99 / 0.997 \times 21.48) \times 0.8937 = \mathbf{0.9240} \text{ Centipoise}$$



Experiment No 8
To determine viscosity of given liquid by Ostwald's Viscometer

- **Result:** Viscosity of given liquid = 0.9240 Centipoises at 25 °C